



Urinalysis reagent for visual or machine reading

# AUTION Sticks

## AUTION Sticks 10V



History of excellence ARKRAY's know-how and expertise has resulted in the development of AUTION Sticks, the most advanced urine test strips.

Urinalysis reagent for visual or machine reading

# AUTION Sticks

AUTION Sticks 10V

## Related products

Compact Urine Analyzer

Pocketchem™ UA

PU-4010

Code: 102678



Versatile Urine Analysis System

AUTION ELEVEN™

AE-4020

Code: 102676



Fully Automated Urine Analyzer

AUTION MAX™

AX-4060



**AUTION Sticks provide both visual and automatic measurement with ARKRAY'S exclusive automated analyzer.**

### Specification

	Measurement principle	Measurement range	Visual measurement time
<b>Glu</b>	Glucose oxidase reaction	50-1000 mg/dl	60 sec.
<b>PRO</b>	Protein-error reaction	15-1000 mg/dl	60 sec.
<b>BIL</b>	Azo-coupling reaction	0.5-6.0 mg/dl	60 sec.
<b>URO</b>	Azo-coupling reaction	2-8 mg/dl	60 sec.
<b>pH</b>	pH indicator	pH 5-9	60 sec.
<b>S.G.</b>	Cation extraction	S.G. 1.000-1.030	60 sec.
<b>BLD</b>	Activity measurement of pseudoperoxidase in hemoglobin	hemoglobin 0.03-1.0mg/dl	60 sec.
<b>KET</b>	Sodium nitroprusside method	5-150 mg/dl	60 sec.
<b>NIT</b>	Griess reaction	nitrite 0.08-0.5 mg/dl	60 sec.
<b>LEU</b>	Activity measurement of esterase in leukocytes	25-500 Leu/μl	90 sec.

**ARKRAY GLOBAL BUSINESS, INC.**

Yousuien-nai, 59 Gansuin-cho,  
Kamigyo-ku, Kyoto 602-0008, JAPAN

TEL 81-75-662-8979 FAX +81-75-431-1202

