



**Abbott**

YOUR FIRST POINT OF CARE

# NYCOCARD™ READER II

NycoCard™ READER II is a multi-assay instrument for point-of-care testing. With NycoCard READER II, healthcare professionals in doctors' offices, clinics and hospitals will obtain accurate, quantitative results within minutes.



# NycoCard™ READER II

## NYCOCARD – AN ESTABLISHED TECHNOLOGY

The NycoCard™ READER II leverages robust technology with proven quality. The instrument is factory calibrated to ensure accurate results. It is a portable, point-of-care instrument designed for rapid and reliable measurements of the NycoCard test. Lot-specific calibration by the operator is not necessary. Only a simple black/white calibration at start-up and the instrument is ready for use.

## BENEFITS AT A GLANCE

- Two tests on one platform: HbA1c and CRP
- Easy to use
- Rapid, quantitative results
- Ideal for measurement of small series
- Local language selection
- Easy maintenance
- Portable, battery-powered (optional)

## SIMPLE TEST PROCEDURE

### 1 Adjusting/ black calibration



### 2 White calibration



### 3 Measure NycoCard™ test



## TECHNICAL DATA

### PRODUCT INFORMATION

Product Code	1116086
Weight	540 g ( Including Reader Pen and Batteries)
Display Type	LCD, 2 × 16 Characters
Communication Interface	RS-232

### DIMENSIONS

Instrument Box	200 mm × 170 mm × 70 mm
Reader Pen	L: 144 mm, D: 29.5 mm

### POWER SUPPLY

Batteries	Rechargeable NiMH Batteries Only, Size AA 1.2 V, Recharge 8–18 hours
Battery Charger	Mascot Type 9725 or Equivalent

### OPERATING CONDITIONS

Temperature	15°C–35°C (The recommended range is 18°C–25°C.)
Relative Humidity	0%–90%. Avoid direct sunlight/too dim light.

## TESTS AND CONTROLS

NYCOCARD™ HbA1c  
NYCOCARD™ CRP

NYCOCARD™ HbA1c CONTROL  
AFINION™ CRP CONTROL

Country Distributor

  
**HEALTH RESOURCES INC.**  
No.13 Sorsogon Street, West Avenue, Quezon City  
Telephone Numbers 87320835 • 87125389  
87314681 • 82570643  
Website : [www.biocare.ph](http://www.biocare.ph)

© 2020 Abbott. All rights reserved. All trademarks referenced are trademarks of either the Abbott group of companies or their respective owners. Any photos displayed are for illustrative purposes only. Any person depicted in such photos is a model. 10001281-07 03/20

  
**Abbott**