



**KAPITOL**

KAPITOL GROUP INTERNATIONAL LTD.

# STATLYTE C PLUS

Electrolyte Analyzer



## Features

- Linux operating system, 800X480 color LCD touch screen.
- High performance electrodes : high precision and wide linearity range, with long shelf life.
- Excellent electrode shielding system, removing external interference and ensuring good stability.
- Real-time detection of each electrode working status, eliminating electrodes drift, ensuring accurate result.
- Possible to connect with a bar-code scanner, accidental power-off protection and automatic fault diagnosis function.
- Support high / low value calibrators and levels of QC. automatic calibration of slope and intercept values. display of whole QC graphs.
- Two sampling modes : automatic sampling with auto sampler or manual sampling.
- Operational software can be upgraded through USB port.
- Automatic liquid level detection and alarms. Real-time diagnostic of system status and reagent pack status information.
- More than 100,000 test results can be stored with convenient sample results review and print-out sample ID.

## Technical Index

Parameter	CV%	Resolution	Range (mmol/L)
K	≤1.0%	0.01mmol/L	0.50~15.00mmol/L
Na	≤1.0%	0.1mmol/L	30.0~200.0mmol/L
Cl	≤1.0%	0.1mmol/L	30.0~200.0mmol/L
Ca <sup>2+</sup>	≤3.0%	0.01mmol/L	0.10~5.00mmol/L
LI	≤5.0%	0.01mmol/L	0.30~3.00mmol/L
pH	≤2.0%	0.01	4.0~9.50 (Unit)
CO <sup>2</sup>	≤3.0%	0.1mmol/L	6.0~50.0mmol/L



No.154, Sec.1, Xinnan Rd., Luzhu District, Taoyuan City, 33858, Taiwan.  
Tel +886-3-3227755 Fax +886-3-3214953 E-mail sales@kapitol-group.com

Country Distributor



No. 13 Sorsogon Street, West Avenue, Quezon City • Telephone (632) 8-7320835 • 8-7314681 • 8-2570643