

LABiTec System Reagents

Reagents	Description	Packing	Test per Box
Prothrombin Time PT	LABiTec -PT- (ISI ~ 1.0) Recombinant thromboplastin; liquid	10 x 2mL 10 x 4mL	200 400
Activated Partial Thromboplastin Time APTT	LABiTec -APTT- Soy phospholipids, ellagic acid, liquid	10 x 2mL 10 x 4mL	400 800
Calcium Chloride CaCl ₂	LABiTec -Calcium Chloride (0.02 M)- for APTT determination, liquid	10 x 4mL	800
Fibrinogen	LABiTec -Fibrinogen Kit- for quantitative determination of Fibrinogen, incl. - LABiTec Thrombin Reagent - LABiTec Imidazole Buffer - LABiTec Kaolin Suspension	5 x 2mL 4 x 25mL 1 x 11mL	200

Controls	Description	Packing	Test per Box
Normal Control	LABiTec -Control Plasma 1- Normal Coagulation Time	10 x 1mL	200
Abnormal Control	LABiTec -Control Plasma 2- Abnormal Coagulation Time	10 x 1mL	200

Calibrators	Description	Packing	Test per Box
Standard Plasma	LABiTec -Standard Plasma- Calibration Plasma for Fibrinogen & PT	10 x 1mL	--

For further details in coagulation and platelet aggregation visit us at www.labitec.de

LABiTec Labor BioMedical Technologies GmbH,
An der Strusbek 6, 22926 Ahrensburg, Germany,
☎ +49-(0)4102-47950, 📠 +49-(0)4102-479535,
✉ info@labitec.de; www.labitec.de



AUTHORIZED DISTRIBUTOR:



BIOCARE HEALTH RESOURCES
64 WEBJET BLDG. COR BMA AVE. QUEZON AVE,
QUEZON CITY, PHILIPPINES
www.biocare.ph